

Infra-Heat Massage Magnetic Cushion




User's Manual

Infra-Heat Massage Magnetic Cushion by Wagan Tech®



User's Manual—Read before using this equipment



Infra-Heat Massage Magnetic Cushion by Wagan Tech®

Thank you for purchasing the Infra-Heat Massage Magnetic Cushion by Wagan Tech®. With normal care and proper treatment it will provide years of reliable service. Please read all operating instructions carefully before use.

WARNING—READ ALL INSTRUCTIONS BEFORE USING

- Always unplug the appliance from the electrical outlet or vehicle immediately after using and before cleaning.
- **DO NOT** reach for an appliance that has fallen into water. Unplug it immediately.
- Keep dry—**DO NOT** operate in a wet, moist environment, or around water.
- **NEVER** use pins or other metallic fasteners with this appliance.
- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use.
- Close supervision is necessary when this appliance used by, on or near children or disabled persons.
- Use this appliance only for its intended use as described in this manual.
- **NEVER** operate this appliance if it has a damaged cord or plug.
- Keep cord away from heated surfaces.
- **NEVER** use while sleeping or fall asleep while using the product.
- **NEVER** drop or insert any object into any opening.
- **DO NOT** operate where aerosol (spray) products are being used or where oxygen is being administered.
- **DO NOT** operate under a blanket or pillow. Excessive heating can occur and cause fire, electrocution or injury to persons.
- This appliance is designed for personal, non-professional use only.

- **DO NOT** use on an infant, child, disabled person, sleeping or unconscious person.
- **DO NOT** use on insensitive skin or on a person with poor blood circulation.
- **DO NOT** crush.
- If you have any concerns regarding your health, consult your doctor before using this product.
- Individuals with pacemakers should consult a doctor before use.
- **DO NOT** use this product for more than 20 minutes at a time.
- Extensive use could lead to the product's excessive heating and shorten the product's life.
- If excessive heating should occur, discontinue use and allow the unit to cool before operating.
- **NEVER** use this product directly on swollen or inflamed areas or skin lesions.
- **DO NOT** use this product as a substitute for medical attention.
- This product should **NEVER** be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls or who has sensory deficiencies in the lower half of their body.

OPERATION

For Vehicle Use:

- Remove the supplied 12V adapter from the package.
- Plug in the output side of the 12V cable into the Infra-Heat Magnetic Cushion.
- Plug the 12V tip into the 12V socket in your vehicle.
NOTE: Be sure to unplug the unit when the vehicle is off to avoid possibly depleting your vehicle's battery.

For Home Use:

- Remove the supplied 110V AC adapter from the package.
- Plug in the output side of the 12V cable into the Infra-Heat Magnetic Cushion.
- Plug the AC adapter into your wall outlet.

To turn on, press the button to the "I" position.

To turn off, press the button to the "O" position

USES

5 Massaging Motors:

- Massages lower and upper back muscles.
- Energizes muscle tissues.
- Loosens tense upper shoulder muscles.
- Perfect for upper thigh and leg muscles.

Heat Wave Therapy:

- Soothes and relaxes tense muscles.
- Increases blood circulation.
- Allows massage to penetrate deeper into muscle tissues.
- Keeps you warm in cold weather

2 Infrared Heat Lights:

- Allows heat to penetrate deeply into muscle tissues to relax and soothe tired muscles.
- Soothing Infrared Heating Points focus in on spots most prone to back tension.

Bio-Magnetic Therapy:

- 10 embedded Bio-Magnets, 950 Gauss each.

Infra-Heat Massage Magnetic Cushion by Wagan Tech®

- Magnetic therapy has been used for centuries by the Chinese.
- Bio-Magnets emit a magnetic field believed to promote blood and oxygen circulation in muscle tissues and acupressure points.

Lumbar Support:

- Stress reducing lumbar support.
- Protects lower spinal discs.
- Promotes better spine alignment.

MAINTENANCE

To Clean:

- Before cleaning, remove the power cord from the outlet and allow unit to cool.
- Wipe with a soft cloth. **DO NOT** use abrasive cleaners.
- Never submerge the unit in any liquid. Keep away from all solvents and harsh detergents.

Storage:

Place massager in a safe, dry and cool place. Avoid contact with sharp edges or pointed objects which might cut or puncture the surface. To avoid breakage, **DO NOT** wrap the power cord around the unit.

WAGAN Corp. Limited Warranty

The WAGAN Corp warranty is limited to products sold only in the United States.

All Wagan Tech products must be registered within (30) days of purchase to activate its warranty. To register your product on our online, please visit <http://tinyurl.com/waganwarranty>. To register your product through mail, please mail in the warranty registration form and a copy of your sales receipt to the address listed on the form. Be sure to keep the original receipt as it will be required when returning a product under the warranty.

Warranty Duration: This product is warranted to the original purchaser for a period of one (1) Year from the original purchase date, to be free of defects in material and workmanship. WAGAN Corporation disclaims any liability for consequential damages. In no event will WAGAN Corporation be responsible for any amount of damages beyond the amount paid for the product at retail.

Warranty Performance: During the above one (1) Year warranty period, a product with a defect will be replaced with a comparable model when the product is returned to WAGAN Corporation with an original store receipt. The replacement product will be in warranty for the balance of the original one (1) Year warranty period.

To return a defective item, please contact WAGAN Corporation at (800) 231-5806 to obtain a Returned Merchandise Authorization number (RMA#), and return instructions. Each item returned will require a separate RMA#. After you have received the RMA# and the return instructions from WAGAN Corporation, please follow the instructions and send the item with PREPAID SHIPPING, along with all of the required documentation, a complete explanation of the problem, your name, address and daytime phone number. WAGAN Corporation will, at its option, replace or repair the defective part.

A Returned Merchandise Authorization number (RMA#) is REQUIRED when sending in any defective item. WAGAN Corporation is not responsible for any item(s) returned without an official Returned Merchandise Authorization number. The item(s) must be returned with prepaid shipping. WAGAN Corporation is not responsible for any shipping charges incurred in returning the item(s) back to the company for repair or replacement. This warranty is void if the product has been damaged by accident, in shipment, unreasonable use, misuse, neglect, improper service, commercial use, repairs by unauthorized personnel or other causes not arising out of defects in materials or workmanship. This warranty does not extend to any units which have been used in violation of written instructions furnished.

Warranty Disclaimers: This warranty is in lieu of all warranties expressed or implied and no representative or person is authorized to assume any other liability in connection with the sale of our products. There shall be no claims for defects or failure of performance or product failure under any theory of tort, contract or commercial law including, but not limited to negligence, gross negligence, strict liability, breach of warranty and breach of contract.



www.wagan.com
customerservice@wagan.com
US & Canada Toll Free: 1.800.231.5806
Tel: +1.510.471.9221

31088 San Clemente Street
Hayward, CA 94544
U.S.A.